



NIH AIDS Reagent Program

20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent: HIV-1 clone ZA012.29

Catalog Number: 11707

Lot Number: 080220

Release Category: C

Provided: 10 μ L purified plasmid DNA, at 1 μ g/ μ L.

Cloning Vector: CMV/R-MCS

Cloning Site: The HIV-1 env/rev cassette was cloned directly into the cloning site of CMV/R-MCS expression vector, in the correct orientation with the CMV promoter. The size of the insert is 2559 bp.

GenBank: EU855133

Description: A PCR fragment containing full-length env and rev genes was derived from the genomic DNA of infected PBMC. The original virus obtained from the AIDS Reagent Program. The env/rev cassette was cloned into CMV/R-MCS expression vector by a directional cloning approach. A single transformed kanamycin resistant E.coli colong was selected and expanded. Recombinant plasmid carries resistance genes for kanamycin.

Special Characteristics: This clone expresses a functional env/rev cassette and can be used to generate psuedotyped infectious virions.

Recommended Storage: -20°C or lower

Contributor: Dr. John Mascola

References: Nature Medicine, volume 13, no. 9, Sept 2007.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

NOTE:

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 clone ZA012.29 (cat#11707) from Dr. John Mascola." Also include the references cited above in any publications.

Scientists at for-profit institutions or who intend commercial use of this reagent must contact Dr. John Mascola, Deputy Director, Vaccine Research Center, National Institute of Health, 40 Convent Drive, Room 4500, Bethesda, MD 20892, Email: jm557m@nih.gov, Tel:(301)594-8487.

Last Updated:

November 07, 2016

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.