

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Anti-HIV-1 RT Monoclonal (MAb 21) Reagent:

Catalog Number: 3483

130167 Lot Number:

С **Release Category:**

Provided: 250 µl ascites fluid.

Host: Balb/c splenocyte x NS-1 myeloma.

Special

Raised against purified HIV-1 reverse transcriptase. Recognizes a central region of the protein, approximately between amino acids 290 and 310. Use at 1:200 for Western Characteristics:

Blot and traditional colorimetric detection and 1:500 and up for chemiluminescent

detection.

Recommended

Storage:

Keep at 4°C for short term storage and -80°C for long term storage. Avoid

freeze-thaw cycles as reagent degradation may result.

Dr. Stephen Hughes. **Contributor:**

Isotype: IgG_1 .

Ferris AL, Hizi A, Showalter S, Pichuantes S, Babe L, Craik CS, Hughes SH. References:

Immunologic and proteolytic analysis of HIV-1 reverse transcriptase structure.

Virology 175:456-464, 1990.

Hizi A, McGill C, Hughes SH. Expression of soluble, enzymatically active, human immunodeficiency virus reverse transcriptase in Escherichia coli and analysis of

mutants. Proc Natl Acad Sci USA 85:1218-1222, 1988.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 11/05/2018 Page 1 of 2

NOTE:	Acknowledgment for	publications should read	"The following reage:	nt was obtained
	/ teltilowicaginent for	publications should read	The following reager	it was obtained

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Anti-HIV-1 RT Monoclonal (MAb 21) from Dr. Stephen Hughes." Also include the references cited above in any publications.

Last Updated November 05, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

Page 2 of 2 REV: 11/05/2018