

## NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

## **DATA SHEET**

Reagent:	HIV-1 Consensus B Gag Peptide Set
Catalog Number:	8117
Lot Number:	130089
Provided:	123 vials, 1.0 mg each, lyophilized.
Description:	All peptides in this set are pooled together and offered as $\#12425$ HIV-1 Consensus B Gag Peptides - Complete Set Pooled.
Sequence:	See Table.
Purity:	>80% as determined by HPLC.
Solubility:	See attached table. Peptides that are difficult to solubilize can almost always be dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to precipitate out of solution.
Solubility:	dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to
Solubility: Special Characteristics:	dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to precipitate out of solution. Please note: Cys and Met containing peptide sequences are highly unstable when stored in DMSO for prolonged periods of time; peptide decay is often observed within days to weeks. If DMSO is required, we recommend freshly-prepared working solutions from freeze dried aliquots of these
Special	dissolved in DMSO. Once a peptide is in solution, the DMSO can be slowly diluted with aqueous medium. Care must be taken to ensure that the peptide does not begin to precipitate out of solution. Please note: Cys and Met containing peptide sequences are highly unstable when stored in DMSO for prolonged periods of time; peptide decay is often observed within days to weeks. If DMSO is required, we recommend freshly-prepared working solutions from freeze dried aliquots of these peptides, and prepared pools.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS. Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 Con B Gag Peptide Set.

Last Updated:

September 10, 2015

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.