

# **Certificate of Analysis for NR-52228**

### Acinetobacter baumannii, Strain MRSN 32076

## Catalog No. NR-52228

This reagent is the tangible property of the U.S. Government.

### **Product Description:**

Acinetobacter baumannii (A. baumannii), strain MRSN 32076 was isolated in 2006 from a human blood sample in Europe as part of a global surveillance program. A. baumannii, strain MRSN 32076 was deposited as sensitive to amikacin, ampicillin/sulbactam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, colistin, gentamicin, imipenem, levofloxacin, meropenem, tetracycline, trimethoprim/sulfamethoxazole and tobramycin. NR-52228 lot 70042852 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70042852 Manufacturing Date: 17MAR2021

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: <a href="mailto:Contact@BEIResources.org">Contact@BEIResources.org</a>. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and
		cream (Figure 1)
Growth at 44°C ± 2°C¹	Growth	Growth
1 day in an aerobic atmosphere on Tryptic		
Soy agar	<b>.</b>	A 4 (1)
Motility  DDL IM Motility Took Modition w/TTC Indicator	Report results	Motile
BBL™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere		
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile <sup>2,3</sup>	71. Dadmaniii	71. Dadmanii (00.070)
Amikacin	Sensitive	Sensitive (8 µg/mL)
Ampicillin/sulbactam	Sensitive	Sensitive (≤ 2 µg/mL)
Cefepime	Sensitive	Sensitive (≤ 2 µg/mL)
Ceftriaxone	Sensitive	Intermediate (16 µg/mL) <sup>4</sup>
Ceftazidime	Sensitive	Sensitive (8 µg/mL)
Ciprofloxacin	Sensitive	Sensitive (≤ 0.25 µg/mL)
Colistin	Sensitive	Sensitive (≤ 0.25 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL)
Imipenem	Sensitive	Sensitive (≤ 1 μg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 1 μg/mL)
Meropenem	Sensitive	Sensitive (≤ 1 μg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 0.5 μg/mL)
Tobramycin	Sensitive	Sensitive (≤ 1 µg/mL)
Tetracycline	Sensitive	Sensitive (4 µg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	99.9% sequence identity to
(~ 1410 base pairs)	A. baumannii, strain MRSN 32076	A. baumannii, strain MRSN 32076
	(GenBank: VHFM01000076.1)	(GenBank: VHFM01000076.1)
Purity	Growth consistent with expected colony	Growth consistent with expected colony
7 days at 37°C in an aerobic atmosphere with	morphology	morphology
and without 5% CO <sub>2</sub> on Tryptic Soy agar		

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898



## **Certificate of Analysis for NR-52228**

TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

<sup>&</sup>lt;sup>1</sup>Growth at 44°C differentiates A. baumannii from A. calcoaceticus and A. pittii, which do not grow at 44°C.

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch 23 JUL 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of VITEK®2 GN82, Sensititre™ GNX2F AST and E-test strips.

<sup>&</sup>lt;sup>4</sup>Susceptibilty result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.