

## **Certificate of Analysis for MRA-201G**

#### Genomic DNA from Plasmodium falciparum, Strain D10

#### Catalog No. MRA-201G

### **Product Description:**


Genomic DNA was extracted from a preparation of *Plasmodium falciparum* (*P. falciparum*), strain D10. MRA-201G lot 63834912 was extracted from a cell culture of BEI Resources MRA-201 seed lot using proprietary technology and vialed in TE buffer (10 mM Tris-HCl and 0.5 mM EDTA, pH 9). Testing was completed on bulk material prior to freezing.

Lot: 63834912 Manufacturing Date: 23NOV2015

TEST	SPECIFICATIONS	RESULTS
Sequencing of Merozoite Surface Protein 2 (MSP2) gene (~ 740 base pairs)	Consistent with <i>P. falciparum</i>	Consistent with P. falciparum
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by Nanodrop Measurement	Report results	~ 0.5 μg in 50 μL per vial (10 μg per mL)
Amount per Vial <sup>1</sup>	Report results	~ 0.5 µg
Functional Activity by PCR Amplification MSP2 locus <sup>1</sup>	600 to 900 base pair amplicon	~ 900 base pair amplicon
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.6 to 2.1	2.0
Protozoan Inactivation Human erythrocytes exposed to 10% of total yield of MRA-201G <sup>2,3</sup>	No parasitemia observed	No parasitemia observed
Mycoplasma Contamination  DNA Detection by PCR  Driver and productions for BCR are succitable.	None detected	None detected

<sup>&</sup>lt;sup>1</sup>Primer sequences and conditions for PCR are available upon request.

Figure 1: Agarose Gel Electrophoresis



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: ~ 673 ng of MRA-201G

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>14 days in complete RPMI culture medium at 37°C in sealed flasks outgassed with blood-gas atmosphere (90% N<sub>2</sub>, 5% CO<sub>2</sub>, 5% O<sub>2</sub>). Complete RPMI culture medium was changed and parasitemia checked every 2 to 4 days.

<sup>&</sup>lt;sup>3</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of *Plasmodium* parasites.



# **Certificate of Analysis for MRA-201G**

/Heather Couch/ Heather Couch

12 AUG 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898