

**Genomic RNA from Zika Virus, MR 766**

**Catalog No. NR-50085**

**Product Description:**

Genomic RNA was isolated from a preparation of cell lysate and supernatant from *Chlorocebus* (previously *Cercopithecus aethiops* kidney epithelial cells (Vero; ATCC® CCL-81™) infected with Zika virus (ZIKV), MR 766 (BEI Resources lot 64254618) using QIAamp® Viral RNA Mini Kit (Qiagen® 52906). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

**Lot: 70041163**

**Manufacturing Date: 20DEC2021**

| TEST   | SPECIFICATIONS                                       | RESULTS   |
|--|--|---|
| <b>Genotypic Analysis</b><br>Sequencing of species-specific region (~ 930 nucleotides)   | ≥ 98% identity with ZIKV, MR 766 (GenBank: AF013415) | 100% identity with ZIKV, MR 766 (GenBank: AF013415) |
| <b>Functional Activity by RT-PCR Amplification<sup>1</sup></b><br>NS5 gene   | ~ 1030 base pair amplicon                            | ~ 1030 base pair amplicon                           |
| <b>Estimated Concentration (post-dilution) by Qubit® Measurement (Viral, Cellular and Carrier)<sup>2</sup></b>                                 | Report results                                       | 12.3 ng per 100 µL (0.12 µg/mL)                     |
| <b>Estimated Amount per Vial<sup>3</sup></b>   | Report results                                       | 12.3 ng   |
| <b>Virus Inactivation</b><br>10% of total yield inoculated on Vero cells and evaluated for cytopathic effect after serial passage <sup>3</sup> | No viable virus detected                             | No viable virus detected                            |

<sup>1</sup>Amplified using iTaq™ Universal SYBR Green One-step Kit (Bio-Rad® 172-5151) with 5 µL of NR-50434 in a 50 µL reaction

<sup>2</sup>Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method

<sup>3</sup>Use of the QIAamp® Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of ZIKV as shown by the absence of cytopathic effect (CPE) after plating the entire extract on virus-susceptible cells for two passages.

/Sonia Bjorum Brower/

Sonia Bjorum Brower

05 APR 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

