

# **Product Information Sheet for HM-115**

## Neisseria flavescens, Strain SK114

# Catalog No. HM-115

## For research use only. Not for use in humans.

#### Contributor:

Guillermo I. Perez-Perez, D.Sc., Associate Professor of Medicine, Departments of Medicine and Microbiology, School of Medicine, New York University, New York, New York, USA

#### Manufacturer:

**BEI Resources** 

## **Product Description:**

Bacteria Classification: Neisseriaceae, Neisseria

Species: Neisseria flavescens

Strain: SK114

<u>Original Source</u>: *Neisseria flavescens (N. flavescens*), strain SK114 was isolated from normal skin of the right arm of a 47-year-old woman.<sup>1,2</sup>

<u>Comments</u>: N. flavescens, strain SK114 (<u>HMP ID 600</u>) is a reference genome for <u>The Human Microbiome Project</u> (HMP). HMP is an initiative to identify and characterize human microbial flora. The complete genome of N. flavescens, strain SK114 was sequenced at the <u>J. Craig Venter Institute</u> (GenBank: <u>ACQV000000000</u>).

Note: HMP material is taxonomically classified by the depositor. Quality control of these materials is only performed to demonstrate that the material distributed by BEI Resources is identical to the deposited material.

*N. flavescens* is an aerobic, Gram-negative diplococcus that is a common inhabitant of the normal human mouth and skin. It is an opportunistic pathogen that has been implicated in rare cases of meningitis, endocarditis and septicemia.<sup>3,4,5,6,7,8</sup>

#### **Material Provided:**

Each vial contains approximately 0.5 mL of bacterial culture in *Haemophilus* Test medium supplemented with 10% glycerol.

<u>Note</u>: If homogeneity is required for your intended use, please purify prior to initiating work.

### Packaging/Storage:

HM-115 was packaged aseptically in cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

## **Growth Conditions:**

Media:

Haemophilus Test medium or equivalent Chocolate agar or equivalent

Incubation:

Temperature: 37°C

Atmosphere: Aerobic with 5% CO<sub>2</sub>

#### Propagation:

- 1. Keep the vial frozen until ready for use, then thaw.
- Transfer the entire thawed aliquot into a single tube of broth.
- 3. Use several drops of the suspension to inoculate an agar slant and/or plate.
- 4. Incubate the tube, slant and/or plate at 37°C for 1 day.

## Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH as part of the Human Microbiome Project: *Neisseria flavescens*, Strain SK114. HM-115."

## Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories (BMBL). Current Edition. Washington, DC: U.S. Government Printing Office.

#### **Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

## **Use Restrictions:**

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

**BEI Resources** 

www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



# **Product Information Sheet for HM-115**

#### References:

- 1. Perez-Perez, G. I., Personal Communication.
- 2. HMP ID 600 (Neisseria flavescens, strain SK114)
- Branham, A. "A New Meningococcus-like Organism (Neisseria flavescens n. sp.) from Epidemic Meningitis." Ann. Intern. Med. 4 (1930): 101.
- Coovadia, Y. M. "Neisseria flavescens Septicaemia with Meningitis. A Case Report." S. Afr. Med. J. 66 (1984): 308-309. PubMed: 6474295.
- Knapp, J. S. "Historical Perspectives and Identification of Neisseria and Related Species." <u>Clin. Microbiol. Rev.</u> 1 (1988): 415-431. PubMed: 3069201.
- Prentice, A. W. "Neisseria flavescens as a Cause of Meningitis." <u>Lancet</u> 272 (1957): 613-614. PubMed: 13407080.
- Sinave, C. P. and K. R. Ratzan. "Infective Endocarditis Caused by *Neisseria flavescens*." <u>Am. J. Med.</u> 82 (1987): 163-164. PubMed: 3799678.
  Wertlake, P. T. and T. W. Williams, Jr. "Septicaemia
- Wertlake, P. T. and T. W. Williams, Jr. "Septicaemia Caused by *Neisseria flavescens*." <u>J. Clin. Pathol.</u> 21 (1968): 437-439. PubMed: 4972296.

 $\ensuremath{\mathsf{ATCC}}^{\otimes}$  is a trademark of the American Type Culture Collection.



BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898