



NIH AIDS Reagent Program

20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent: pV16970_2

Catalog Number: 11809

Lot Number: 110247

Release Category: A

Provided: 5 µg of dried purified DNA stabilized in DNastable *PLUS*

Cloning Vector: See pNL4.3 GenBank accession AF324493.1 for the vector sequence. The vector contains an AmpR marker.

Cloning Site: Restriction maps match those of pNL4.3 with the exception of a novel Bal I site at nucleotide position 4552 and PflM 1 site at nucleotide position 5302 that were removed and a PflM 1 site added at nucleotide position 3491. The unique Apa I site upstream of protease at nucleotide position 2009 is used as a 5' restriction site; double-digestion with Apa I/Bal I generates a protease-deleted linearized vector. Inserts are amplified from clinically-derived viral cDNA using Pfu; Purified PCR products are digested and ligated into the vector using cycling T4 ligation, and the ligation product is grown in Stb13 E. coli in LB+carbenicillin, then plated and sequenced. The cloned region contains the entire protease open reading frame; the flanking regions cloned along with protease include 3' end of gag with the gag cleavage site, and the 5' end of RT.

GenBank: [KC109801](#)

Host Strain: Plasmids can be propagated in STBL2 cells and grown at 37°C. Larger plasmids may benefit from growth at 30°C. This construct may also be grown in other competent cells.

Description: This clone (clone V16970_2, PtID 1329) is part of a panel containing each of the canonical protease inhibitor-resistance pathways in clinically derived HIV-1 clones. In contrast to site-directed mutants, the mutations are present in their naturally occurring genetic contexts, which may include known accessory drug-resistance mutations, as well as changes at positions that are not currently known to be associated with drug resistance. As the clones are also infectious and replication-competent, they can be used for in vitro

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

susceptibility testing or new protease inhibitors (PIs). PIs that are active against these clones are likely to retain activity against the most clinically relevant, or possibly all, protease inhibitor-resistant variants. Researchers can also create their own recombinant viruses using the pNLFPB vector. A protocol is available upon request (please contact Shafer Lab)

Please see this [TABLE](#) for the other clones in this panel.

Special Characteristics:

This construct is 14825 bp including the insert.

[Sequence file lot 110247](#)

This reagent is currently being provided as dried purified DNA stabilized in DNastable *PLUS*. Please see the notice for additional information and the protocol for reconstitution of dried DNA reagents. [Dried DNA Notice](#)

Recommended Storage:

Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier bag.

Contributor:

Robert Shafer

References:

Prototypical Recombinant Multi-Protease-Inhibitor-Resistant Infectious Molecular Clones of Human Immunodeficiency Virus Type 1 Vici Varghese, Yumi Mitsuya, W. Jeffrey Fessel, Tommy F. Liu, George L. Melikian, a David A. Katzenstein, Celia A. Schiffer, Susan P. Holmes, and Robert W. Shafer. *Antimicrob Agents Chemother.* 2013 September; 57(9): 4290-4299. doi: 10.1128/AAC.00614-13. [PMID: 23796938](#).

NOTE:

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH from Dr Robert Shafer: Cat# 11809, pV16970_2" Also include the reference cited above in any publications.

Last Updated:

January 29, 2019

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.