

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Reagent: ★ HIV-1 91/US/2 (GS 005) Virus

Catalog Number: 7687

098297 Lot Number:

Release Category: D

Provided: 1 mL cell-free virus

 $TCID_{50} = 6.3 \times 10^2 / mL$

p24 = 61 ng/mL

Host Strain: Human PBMC

Sterility: Negative for mycoplasma, bacteria and fungi

Description: Subtype B

Special

Characteristics:

Part of a a large panel of well-characterized, cultured HIV-1 standard isolates

representing subtypes A through CRF02_AG.

Click here to see the table of reagents that make up this panel

Genbank number: L14574

Recommended

Storage:

Liquid nitrogen

Contributor: Dr. Nelson Michael

References: Jagodzinski LL, Wiggins DL, McManis JL, Emery S, Overbaugh J, Robb M, Bodrug S,

Michael NL. Use of calibrated viral load standards for group M subtypes of human immunodeficiency virus type 1 to assess the performance of viral RNA quantitation tests.

J Clin Microbiol 38:1247-1249, 2000.

Michael NL, Herman SA, Kwok S, Dreyer K, Wang J, Christopherson C, Spadoro JP,

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 08/10/2018 Page 1 of 2 Young KK, Poionis V, MICCUTCHAIN FE, Carr J, Mascola JK, Jagodzinski LL, Kodd MIL. Development of calibrated viral load standards for group M subtypes of human immunodeficiency virus type 1 and performance of an improved AMPLICOR HIV-1 MONITOR test with isolates of diverse subtypes. *J Clin Microbiol* **37**:2557-2563, 1999.

Vahey M, Nau ME, Barrick S, Cooley JD, Sawyer R, Sleeker AA, Vickerman P, Bloor S, Larder B, Michael NL, Wegner SA. Performance of the Affymetrix GeneChip HIV PRT 440 platform for antiretroviral drug resistance genotyping of human immunodeficiency virus type 1 clades and viral isolates with length polymorphisms. *J Clin Microbiol* **37**:2533-2537, 1999.

NOTE: For commercial use, please contact: Dr. Nelson L. Michael, Director, Division of

Retrovirology, Walter Reed Army Institute of Research. Email:

nmichael@hivresearch.org, Tel: 301-792-3430, http://www.hivresearch.org.

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 91/US/2 (GS 005) Virus from Dr. Nelson Michael." Also include the references cited above in any

publications.

Last Updated: August 10, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 08/10/2018 Page 2 of 2