

Product Information Sheet for NR-56637

Klebsiella pneumoniae, Strain BK34907

Catalog No. NR-56637

For research use only. Not for use in humans.

Contributor and Manufacturer:

ATCC®

Product Description:

Bacteria Classification: Enterobacteriaceae, Klebsiella

Species: Klebsiella pneumoniae

Strain: BK34907

Original Source: Klebsiella pneumoniae (K. pneumoniae), strain BK34907 was isolated in 2011 from a urine sample in

the United States.1

Comments: K. pneumoniae, strain BK34907 was deposited as part of the Global Priority Superbugs Collection. NR-56637 was deposited as resistant to amoxicillin/clavulanate, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, doripenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and trimethoprim/sulfamethoxazole.

K. pneumoniae is a Gram-negative enterobacterium that is a major cause of nosocomial infections of the urinary and respiratory tracts. Due to the extensive spread of antibiotic-resistant strains, especially extended-spectrum β -lactamase (ESBL)-producing strains, there has been renewed interest in Klebsiella infections. 2,3,4

Material Provided:

Each vial contains approximately 0.3 mL of bacterial culture in Tryptic Soy broth supplemented with 10% glycerol.

<u>Note</u>: If homogeneity is required for your intended use, please purify prior to initiating work.

Packaging/Storage:

NR-56637 was packaged aseptically in cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

Media:

Nutrient broth or Tryptic Soy broth or equivalent

Nutrient agar or Tryptic Soy agar or Tryptic Soy agar with 5% defibrinated sheep blood or equivalent

Incubation:

Temperature: 37°C Atmosphere: Aerobic

Propagation:

- 1. Keep vial frozen until ready for use, then thaw.
- Transfer the entire thawed aliquot into a single tube of broth.
- Use several drops of the suspension to inoculate an agar slant and/or plate.
- 4. Incubate the tube, slant and/or plate at 37°C for 1 day.

Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Klebsiella pneumoniae*, Strain BK34907, NR-56637."

Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories (BMBL). Current Edition. Washington, DC: U.S. Government Printing Office.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale. This material may be subject to third party patent rights.

References:

- 1. McGann, P., Personal Communication.
- Lascols, C., et al. "Increasing Prevalence and Dissemination of NDM-1 Metallo-β-Lactamase in India: Data from the SMART Study (2009)." <u>J. Antimicrob.</u> <u>Chemother.</u> 66 (2011): 1992-1997. PubMed: 21676902.

BEI Resources

www.beiresources.org

E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898



Product Information Sheet for NR-56637

- Ramirez, M. S., et al. "Multidrug-Resistant (MDR) Klebsiella pneumoniae Clinical Isolates: A Zone of High Heterogeneity (HHZ) as a Tool for Epidemiological Studies." Clin. Microbiol. Infect. 18 (2012): E254-E258. PubMed: 22551038.
- Podschun, R. and U. Ullmann. "Klebsiella spp. as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors." Clin. Microbiol. Rev. 11 (1998): 589-603. PubMed: 9767057.

ATCC[®] is a trademark of the American Type Culture Collection.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898